



Authorized Distributor



Increasing well-being through digital dentistry.

Unlock the Full Potential of Your Smile Journey

LuxCreo Rejuv™ is your secret weapon for a successful treatment.

Welcome to Your 4D Aligner™ Journey!

Q Is my 4D Aligner™ an oral hygiene indicator?

Yes. Discoloration or odor is a warning sign.

It signals plaque and bacteria buildup, increasing your risk of cavities.



Q Why do my aligners need a "recharge"?

Daily wear can weaken your aligner's shape and force.

The Rejuv™ restores this power to keep your treatment on track, while reversing any discoloration or odor that may appear during use.



Q So, what does the Rejuv™ actually do?

Restores Power: In 10 min, a warm bath restores the aligner's effective shape and force.

Deep Cleans: Tablets kill 99.9% of bacteria that brushing misses.



Is it Easy to Use?

Yes! It's a simple process.

► Rejuvenate Power

In 3 Simple Steps: Integrating the Rejuv™ into your routine is effortless.



► Master Hygiene

Complete Your Daily Oral Care Toolkit:

For the best results, we strongly recommend combining the Rejuv™ with these daily habits.



Mouthwash



Electric Brush



Regular Flossing



Water Pick

Prevent Problems Before

They Start: Twice a day rinse with **OXYGEN-POWERED** mouthwash.



TheraBreath



Colgate Peroxyl



Heritage Store



Authorized Distributor



Maximize Your Investment in a New Smile

Ask your dental team how the LuxCreo Rejuv™ and a great daily routine lead to faster, safer, and more hygienic results.

+1 (888) 852-1755

info@3dnadental.com

www.3DNADental.com

Rejuv™ FAQ

