

COMPATIBILITY

Pac-Dent is dedicated to expanding this list of compatible devices as we believe in giving you the freedom of choice in using your preferred devices with all of our resins.

Compatible Printers	Rodin [™] Sculpture Ceramic Nanohybrid	Rodin [™] Denture Base	Rodin [™] Clear Bite Splint	Rodin [™] Custom Bite Tray	Rodin [™] Surgical Guide	Rodin [™] Try-In Denture	Rodin [™] Ortho IBT	Rodin [™] Model	Rodin [™] Castable
Ackuretta [®] SOL	✓	✓	✓	✓	✓	✓	✓	✓	⌚
Asiga [®] MAX (385 nm)	✓	✓	✓	✓	✓	✓	⌚	✓	⌚
Asiga [®] Pico 2 (385 nm)	✓	✓	✓	✓	✓	✓	⌚	✓	⌚
Asiga [®] Pro 4K (385 nm)	✓	✓	✓	✓	✓	✓	⌚	⌚	⌚
Phrozen [®] Sonic Mini 4K	✓	✓	●	✓	✓	✓	✓	⌚	⌚
Phrozen [®] Sonic 4K	✓	✓	●	⌚	✓	✓	✓	⌚	⌚
Phrozen [®] Sonic XL 2022	✓	✓	●	⌚	⌚	⌚	⌚	⌚	⌚
Carbon [®] 3D Printers	⌚	⌚	⌚	✓	⌚	⌚	⌚	⌚	⌚
FlashForge [®] Focus 6K	✓	✓	⌚	✓	✓	✓	✓	✓	⌚